




**Be Right™**



## **Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes**

**Product #:** 2497900  
**SGD Price:** Contact Hach  
**Available**

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

---

### **Specifications**

Diameter:	13 mm
Holes:	90
Quantity:	1 each
Test Tube Arrangement:	6 x 15
Volume / package size e Post Text:	each