



Be Right™



## 6-Branch Stainless Steel Manifold

Product #: 2521406  
SGD Price: Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Sartorius Combisart® system enables you to select the optimal hardware and consumables for your needs in microbiological analysis or particle count in quality assurance. Combisart® features a modular design and field-proven standard accessories to make your choice easier. At the heart of the Combisart® system is a stainless steel manifold designed to accommodate Biosart® 100 Monitors and Biosart® 250 Funnels. Used in conjunction with the Stainless Steel Support, cat# 2521410. Autoclavable.

---

### Specifications

Filter Holder:	Removable
Filtration Area:	12.5 cm <sup>2</sup>
Number of Filter Holders:	6
Region:	IN, SG, TH, US
Suitable Membrane Filter Diameter:	47/50 mm