



**Be Right™**



## **Glass Erlenmeyer Flask, 250 mL, with Screw Cap, 6 pack**

**Product #:** 2089676  
**SGD Price:** Contact Hach  
**Ships within 1 week**

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### **Specifications**

**Capacity:** 250 mL  
**Units:** 6/pk