




**Be Right™**



## **Beaker, Griffin, Low Form, Glass, 1000 mL**

<b>Product #:</b>	<b>50053</b>
SGD Price:	Contact Hach
Available	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### **Specifications**

Body material:	Glass
Capacity:	1000 mL
Material:	Glass, Griffin, Low Form
Units:	each