



Be Right™



Soxhlet Thimbles, Cellulose, 22x100mm, 25/PK

Product #:	2966600
SGD Price:	Contact Hach
Ships within 1 week	

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Specifications

Material:	Cellulose
Pore Size:	10 µm
Quantity:	25 /pk