



Be Right™



QL Filters, Pore Size 8 - 12 μ m, Diameter 125 mm, Circular, 100/pk, Qualitative

| | |
|-------------------|--------------|
| Product #: | 50657 |
| SGD Price: | Contact Hach |
| Available | |

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifications

| | |
|-----------|-----------|
| Diameter: | 125 mm |
| Material: | Cellulose |
| Quantity: | 100 /pk |