



Be Right™



Stirrer, Lab, Non-programmable, 6 Paddle, 120V, Jar Tester

Product #:

2703800

 **Hazardous**

SGD Price:
Available

Contact Hach

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Standard model for either four- or six-paddle lab stirrer. Features digital readout, stainless steel paddles placed 6 inches apart, regulated variable speeds of all paddles simultaneously, from 1-300 rpm, powder-coated steel uniframe chassis, built-in illuminator, anti-glare curtain, and dust cover. The unit holds either round 1000 mL mixing beakers or square 2 L Wagner Floc Jars. Beakers and Wagner Floc Jars are sold separately.

Quality craftsmanship

Built on a powder-coated steel uni-frame chassis

Illuminator

Provides reliable, diffused sample illumination

Versatile programming

High accuracy and repeatability