



Be Right™



Manifold, Filter Funnel, Aluminum

| | |
|----------------------|----------------|
| Product #: | 2486100 |
| SGD Price: | Contact Hach |
| Ships within 2 weeks | |

Feature 3-place filtration with independent vacuum control. Handles 25 or 47 mm filter funnels for processing multiple samples simultaneously. Designed for efficient analysis of microbiological or particulate contamination. Aluminum: Anodized aluminum body; chrome-plated brass valves; stainless steel drain plugs. Autoclavable. Dimensions (D x W x H) - 6 x 16 x 6.25" Polyurethane: High-impact polyurethane body; glass-filled polypropylene valves; nylon drain plugs. Dimensions (D x W x H) - 4 x 11 x 6"

Specifications

Quantity: 1 each