



**Be Right™**



## Microfil V Filtration Device

<b>Product #:</b>	<b>2508900</b>
SGD Price:	Contact Hach
Available	

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 72 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45  $\mu\text{m}$  pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

---

### Specifications

Capacity:	100 mL
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 $\mu\text{m}$