



Be Right™



## Hach Modified Colitag Presence-Absence Media, 100/pk Snap-Pack

**Product #:** 8888004  
**SGD Price:** Contact Hach  
**Available**

### Hach Colitag presence-absence tests deliver the fastest EPA approved results in the market

Hach Modified Colitag™ is the fastest presence-absence test in the market. This allows drinking water facilities to begin a test at 5:00 p.m. and get EPA-approved results by 9:00 am the following morning. The product also provides unparalleled flexibility by offering two testing parameters in one product.

#### The Fastest Presence-Absence Test on the Market at 16 Hours

Fastest EPA-approved test for coliforms and *E. coli*. Test as late as 5:00 PM and still have your results by 9:00 AM the next day.

#### One Product Supports 16- or 22-hour Test Methods

Allows for a "standard" test that requires less set-up, or a "rapid" test when necessary. One media for both tests reduces complexity in certification and inventory.

#### Results are Valid for 48 Hours

Saturday's tests can be read on Monday—no need to work all weekend. Results stay valid for 800% longer than the competition.

---

### Specifications

Quantity: 100 /pk

---

### Required Accessories

- Sample Vial, 120 mL with Thiosulfate and Custody Tie, 200/pk (Item 2634001)
- Incubator, Water Bath, 110 Vac (Item 2616300)
- Incubator, Culture, Low-Profile 110 V (Item 2619200)
- Comparator for use with Hach Modified Colitag Presence-Absence Media (Item 8888001)
- Hach Modified Colitag Sterile 120 mL Sample Bottles, 100/pk, Shrink-banded, Polystyrene (Item 8888006)