



**Be Right™**



## **Low Range Orthophosphate Chemkey® Reagents (box of 25)**

**Product #:** 8636600  
**SGD Price:** Contact Hach

**Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.**

Low Range Orthophosphate Chemkeys, range 0.20 - 4.00 mg/L PO<sub>4</sub>. For SL1000 Portable Parallel Analyser (PPA).

**Current SL1000 instrument customers who want to test Orthophosphate need to perform a software update that can be located on the instrument download section.**

---

### **Specifications**

Parameter:	Orthophosphate
Platform :	Chemkeys
Quantity:	pk/25
Range:	0.20 - 4.00 mg/L PO <sub>4</sub>
Shelf Life:	18 months