



Be Right™



High Range Total Alkalinity Chemkey® Reagents (box of 25)

Product #: 8636100
SGD Price: Contact Hach

Use with the SL1000 Portable Parallel Analyser (PPA) for dramatically streamlined water quality testing.

High Range Total Alkalinity Chemkeys, range 200 - 700 mg/L CaCO₃. For SL1000 Portable Parallel Analyser (PPA).

Current SL1000 instrument customers who want to test Alkalinity need to perform a software update that can be located on the instrument download section.

Specifications

Method:	10283
Method Name:	Colorimetric
Parameter:	Alkalinity
Platform :	Chemkeys
Quantity:	pk/25
Range:	200 - 700 mg/L CaCO ₃
Shelf Life:	9 months